Extractor Aspiration Catheter

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REF	Catalogue number	LOT	Batch code
	Do not re-use		Caution
	Do not use if package is damaged		Temperature limit
	Sterilized using ethylene oxide		Keep away from sunlight
	Use by date		Date of manufacture
	Consult instructions for use		Manufacturer
	Keep dry		Do not resterilize

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


MS-P953/IFU-DS1003(Rev.4,2401)

A CAUTION

The Extractor Aspiration catheter procedure should be performed by physicians with adequate training in the use of the device.

B DEVICE DESCRIPTION

The Extractor Aspiration catheter is a dual lumen rapid exchange catheter with related accessories. The smaller wire lumen of the catheter is able to accommodate guidewires that are $\leq 0.014"$ / 0.36mm in diameter. The larger aspiration lumen comes preloaded with a core wire that resists kinking during placement but is removed to allow for the removal of thrombus by aspiration. The catheter have maximum outer diameter of 1.70, 1.90 and 2.10 mm by size and two working length of 140 and 150 cm. The catheter has a radiopaque marker band (1mm) located 3.35mm from the distal tip. The catheter has three (3) non-radiopaque positioning marks located approximately 90cm, 100cm, and 110cm proximal of the distal tip. Catheter has a distal hydrophilic coated section of 40cm. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of the catheter to the included extension line, stopcock and syringes, as well as a self-sealing stylet port. Two 180 micron filter baskets are included for assistance in filtering the blood removed during the procedure for laboratory analysis of any thrombus.

C INDICATIONS

The Extractor Aspiration catheter is indicated for the removal of fresh, soft emboli and thrombi from vessels in the coronary and peripheral vasculature.

D CONTRAINDICATIONS

The Extractor Aspiration catheter is contraindicated for use in patients with:

- Cardiogenic shock
- No identifiable culprit lesion
- Coronary morphology ineligible for thrombectomy (for example, very tortuous vessels or severe calcification)
- Co-morbidity with limited life expectancy of less than 6 months
- Contraindications for treatment heparin, aspirin, or thienopyridines
- Pregnancy

E ⚠ WARNINGS

Do not use after the "Use-by date" specified on the label. Do not use the Extractor Aspiration catheter for the delivery or infusion of diagnostic, embolic or therapeutic materials into blood vessels as it has not been designed for these uses. The core wire shall be removed before initiation of thrombus aspiration. Do not perform high pressure contrast injections around the Extractor Aspiration catheter within the minimum sized guide catheter. High pressure contrast injection may damage the Extractor Aspiration catheter, making it difficult to remove from the guide catheter. There will be insufficient room within the minimum guide catheter sizes to deliver any contrast around the Extractor Aspiration catheter. Do not use the syringes, extension line, stopcock or filter basket inside the human body. If flow into the syringe stops or is restricted, do NOT attempt to flush the extraction lumen while the catheter is still inside the patient's vasculature. Intravascular thrombus delivery, thromboembolic event and/or serious injury or death may result. Remove the catheter and, outside of the patient, either flush the extraction lumen or use a new catheter. Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire tip, damage to the catheter, or vessel perforation. The device is provided sterile for single use only. Reuse of single-use device creates a potential risk of patient or user infections. Contamination of the device may lead to illness or serious patient injury. Do not expose the catheter to organic solvents, e.g. alcohol.

F ⚠ PRECAUTIONS

Do not use the Extractor Aspiration catheter less than 2.0 mm in diameter of blood vessel. Do not use the Extractor Aspiration catheter for the removal of fibrous, adherent or calcified material (e.g. chronic clot, atherosclerotic plaque). Do not use the Extractor Aspiration catheter in the venous system. Do not use the Extractor Aspiration catheter in the cerebral vasculature. Do not use the Extractor Aspiration catheter if the packaging has been damaged. Inspect the catheter prior to use for any

bends or kinks. Do not use a damaged catheter because vessel damage and/or inability to advance or withdraw the catheter may occur. Check that all fittings are secure so that air is not introduced into the extension line or syringe during extraction. Exercise care in handling of the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking. When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response. Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the catheter. The Extractor Aspiration catheter has not been evaluated for use in the venous vessels.

G ADVERSE EFFECTS

The following adverse effects are generally associated with catheterization procedures and may occur when using the Extractor :

- local or systemic infection
- local hematomas
- intimal disruption
- arterial dissection
- perforation and vessel rupture
- arterial thrombosis
- distal embolization of blood clots and plaque
- arterial spasm
- arteriovenous fistula formation
- catheter fracture with tip separation and distal embolization

H CLINICAL PROCEDURE

The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. Each Extractor Aspiration catheter includes the following components:

- Single-use, disposable catheter with core wire
- Extension line (30cm) with attached stopcock
- 2 - syringes (30ml) with locking plunger
- 2 - filter basket (180 micron mesh) for filtering thrombus after extraction

Other materials required but not provided are:

- Guiding catheter with the appropriate sized I.D. to cannulate the vessel
- Rotating hemostatic valve (RHV)
- Guidewire with diameter $\leq 0.014"$ / 0.36mm
- Sterile heparinized saline (for system flushing)

I PREPARATIONS FOR USE

1. Using sterile technique, open the pouch and transfer the tray into the sterile field.
2. Flush the carrier tube completely with heparinized saline to activate the hydrophilic coating.
3. Remove the catheter from the carrier tube and inspect for any bends or kinks. Remove the packaging mandrel from the wire lumen.
4. Thoroughly flush the catheter and wire lumen with heparinized saline. 5. Draw 5ml of heparinized saline into the 30ml syringe. Connect the syringe to the stopcock, and connect the attached extension line to the catheter. Flush the entire connection to remove all air from the catheter, extension line, stopcock and syringe. Turn the stopcock to the "Off" position.
6. With the stopcock in the "Off" position, pull back the plunger on the 30ml syringe to the desired amount of extraction volume. Twist the plunger to lock the syringe in the vacuum position.
Precaution: Check that all fittings are secure so that air is not introduced into the extension line or syringe

J DEPLOYMENT PROCEDURE

Review the schematic drawing of the Extractor Aspiration catheter below for terminology used in the following deployment procedure. Deploy the Extractor™ Aspiration catheter according to the following steps:

1. Cannulate the vessel using the appropriate guidewire and guiding catheter with attached RHV. Flush the guiding catheter and RHV using standard technique.
NOTE: A correctly sized guide catheter must be selected to accommodate free movement of the Extractor™ Aspiration catheter.
2. Back-load the Extractor Aspiration catheter onto the guidewire. Advance the catheter on the guidewire until the wire exits the opening in the wire lumen.

3. Open the RHV thumbscrew and introduce the catheter, being careful to keep the guidewire in the guidewire lumen slot of the catheter. Tighten the Oring valve around the catheter just enough to catheter movement.
4. Continue to advance the catheter over the guidewire to the selected vascular site, using fluoroscopic guidance. Stop advancement of the catheter if any resistance is encountered.

Warning: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

5. After fluoroscopically confirming catheter position (proximal to the thrombus), remove the core wire. Open the stopcock to begin extraction when the hub of catheter is connected to male luer of extension line. Slowly advance the catheter distally away from the guiding catheter. Blood will enter the syringe until all of the vacuum is gone. Should aspiration not begin filling the syringe within 5 seconds remove the catheter without releasing the vacuum. Outside of the patient, either flush the extraction lumen or use a new catheter.

NOTE: During performance testing, the catheter demonstrated the ability to evacuate fluid at a minimum rate of 0.20ml/s.

Warning: If flow into the syringe stops or is restricted, DO NOT attempt to flush the extraction lumen while the catheter is still inside the patient's vasculature. Intravascular thrombus delivery, thromboembolic event and/or serious injury or death may result. Remove the catheter and, outside of the patient, either flush the extraction lumen or use a new catheter.

NOTE: If air is noted in the syringe during extraction, a leak may be present in the system. Turn the stopcock "Off", tighten all luer connections, remove all air from the syringe, and repeat the extraction. If air is still noted, remove the catheter, obtain a new catheter and repeat the procedure.

6. After completing the extraction process, turn the stopcock to the "Off" position and remove the catheter, or attach second syringe and repeat extraction.
7. Blood and thrombus extracted into the syringe may be filtered for subsequent laboratory analysis using the 180 micron filter basket. To use the filter basket, wet the mesh to allow fluid flow, and filter the extracted blood through the filter basket.

K STORAGE

Store at room temperature (1°C-30°C) in dry and dark places.

L LIMITED WARRANTY

GENOSS Co., Ltd. warrants that the Extractor™ Aspiration catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by GENOSS Co., Ltd. to be defective in workmanship or materials. GENOSS Co., Ltd. shall not be liable for any incidental, special, or consequential damages arising from the use of the Extractor™ Aspiration catheter.

Damage to the product through misuse, alteration, improper storage, or improper handling shall void this limited warranty. No employee, agent, or distributor of GENOSS Co., Ltd. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against GENOSS Co., Ltd.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF GENOSS CO., LTD.

M SYMBOLS

REF	Catalogue number	LOT	Batch code
	Do not re-use		Caution
	Do not use if package is damaged		Temperature limit
	Sterilized using ethylene oxide		Keep away from sunlight
	Use by date		Date of manufacture

	Consult instructions for use		Manufacturer
	Authorized representative		Do not sterilize
	Keep dry		

N PRODUCT SPECIFICATION


When using the products, please refer to the sizes which are followed in the table 1. And please confirm firstly the maximum diameter of catheter, compatible guide catheter and the diameter of the guidewire.

[table-1]

Catheter Model	Max. Guidewire Diameter	Min. Guide Catheter I.D.
GAC-6-110	0.014" (0.36mm)	0.070" (1.78mm)
GAC-7-130	0.014" (0.36mm)	0.080" (2.03mm)
GAC-8-150	0.014" (0.36mm)	0.089" (2.26mm)

The sizes of the catheter and core wire A and B please refer to the table 2.

[Catheter]

[Core wire]


[table-2]

Catheter Model	Catheter	Core Wire
	A	B
GAC-6-110	1400mm	1460mm
GAC-7-130	1400mm	1460mm
GAC-8-150	1500mm	1560mm

Accessories, please refer to the following

1. Stopcock Valve: 1EA
2. Extension line: 1EA
3. Pressure Lock Syringe : 2EA
4. Filter basket : 2EA
5. One set of extension line, stopcock valve, pressure lock syringe is packaged with the concluded states. Pressure lock syringe are packaged separately. As shown in the figure-1 below.

[figure-1]

MS-P953/IFU-DS1003(Rev.4,2401)