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REF	Catalogue number	LOT	Batch code
	Do not re-use		Caution
	Do not use if package is damaged		Temperature limit
STERILIZED	Sterilized using ethylene oxide		Keep away from sunlight
	Use by date		Date of manufacture
	Consult instructions for use		Manufacturer
	Do not re-sterilize		Keep dry

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A DESCRIPTION

The GENOSS PTA Balloon Catheter is an Over the Wire (OTW) peripheral balloon catheter, specifically designed for Percutaneous Transluminal Angioplasty in stenosed vessel segments. The balloon is designed to inflate to a controlled diameter and length when inflated. The Balloon Catheter has two platinum iridium radiopaque marker bands. Two radiopaque marker bands are located at both ends of the balloon, to facilitate fluoroscopic visualization and define the length and location of the balloon. The balloon catheter includes a tapered soft-tip to PTA Balloon Catheter are 50cm, 80cm and 130cm. The catheter has a dual lumen shaft. The dual lumen shaft is branched at the proximal end so that one lumen forms the entrance to the central lumen for the guidewire, of diameters 0.035"(0.89mm) or smaller to facilitate advancement of the balloon catheter through the stenosis segment to be dilated, while the other lumen is used to inflate and deflate the dilatation balloon with a mixture of contrast medium and saline solution. The balloon catheter is provided with hydrophilic coating that is activated when wet.

B INDICATIONS FOR USE

A balloon catheter used for percutaneous transluminal angioplasty of peripheral vascular iliac artery, femoral artery, femoral-iliac artery, is used in a closed lesions due to congenital or synthetic arteriovenous dialysis PTA and lower popliteal artery and renal artery.

C CONTRAINDICATIONS

Typically uses a balloon catheter contraindications are as follows
1) Can not pass the lesion along with Guide wire.
2) Do not use in the central circulation system such as coronary and supra-aortic artery.

D WARNINGS

This device is designed and intended for single use only. Do NOT reuse, reprocess or re-sterilize. GENOSS will not be responsible for any direct, incidental or consequential damages resulting from re-sterilization and /or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Do not use if sterile barrier is damaged. If damage is found, call our GENOSS representative. Before operate, use heparinized saline for flushing to prevent clotting. After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy. Prior to procedure, the balloon catheter should be examined to verify functionality and ensure that its size is suitable for the specific procedure for which it is to be used. The inflated diameter of the balloon should never exceed the original diameter of the vessel proximal and distal to the lesion. Balloon pressure should not exceed the rated burst pressure. Use of a pressure monitoring device is mandatory to prevent over-pressurization. Use only an appropriate balloon inflation medium [e.g. 50:50 mixture by volume of contrast medium and saline]. Never use air or any gaseous medium to inflate the balloon. Do not expose the balloon catheter to organic solvents, e.g., alcohol, ionizing radiation or ultraviolet light. Rotate inventory so that catheters are used prior to the expiration date on package label. When the catheter is in the body, it should be manipulated while under sufficient and/or high-quality fluoroscopy. Do not advance the catheter unless the balloon is fully deflated under vacuum. If strong resistance is met during manipulation, stop the procedure and determine the cause of the resistance before proceeding. The short-term and long-term biological effects of balloon diameters larger than the original vessel diameters are not known.

E PRECAUTIONS

Do not use the balloon catheter if the outer or the inner package is damaged or opened. Use the balloon catheter before the "shelf-life" date specified on the package. Before angioplasty, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used. Use only guide wires with a 0.035"(0.89mm) maximum diameter. Exercise care during handling to reduce the possibility of accidental breakage, bending or kinking of the balloon catheter shaft. If strong resistance is met during manipulation, stop the procedure and determine the cause of the resistance before proceeding. Only physicians thoroughly trained and educated in the performance of Percutaneous Transluminal Angioplasty (PTA) should use the balloon catheter. PTA should only be performed at hospitals where emergency peripheral artery bypass graft surgery is readily available.

F POTENTIAL COMPLICATIONS/ ADVERSE EFFECTS

Possible complications include, but are not limited to:

- Abrupt closure
- Access site hematoma
- Aneurysm
- Angina
- Arrhythmias
- Arteriovenous fistula
- Bleeding complications, which may require transfusion
- Cerebral ischemia / transient ischemic attack (TIA)
- Death
- Embolism (air, tissue, thrombotic, systemic or device component)
- Fever / pyrogenic reaction
- Hypersensitivity or allergic reaction to contrast agents and drug reactions
- Hypertension / hypotension
- Infection
- Ischemia, including tissue ischemia, steal syndrome, and necrosis
- Leg edema
- Myocardial ischemia or infarction
- Nausea and vomiting
- Neuropathies or nerve injury
- Occlusion
- Organ failure (single, multiple)
- Pain
- Palpitations
- Pseudoaneurysm
- Renal failure / insufficiency
- Restenosis
- Stroke / cerebrovascular accident (CVA)
- Vascular complications, including entry site, which may require vessel repair
- Vascular thrombosis
- Vessel injury, e.g. dissection, perforation
- Vessel spasm

G PREPARATION TECHNIQUE

Prior to use, carefully examine the unit to verify that the catheter or sterile package has not been damaged in shipment. After moistening the balloon protection sheath with physiological saline solution, carefully remove the balloon protection sheath not to damage the balloon. In case of resistance, twist the protection sheath with one hand, holding the shaft with the other.

CAUTION: Do not use with Lipiodol contrast media, or other such contrast media, which incorporate the components of these agents.

TO DISPLACE AIR:

1. Flush out the lumen for the guidewire by connecting a syringe dilled with saline solution. Start flushing until liquid rinses out of guidewire lumen at the distal tip of the balloon catheter.
2. Connect a syringe containing contrast medium to the Luer connector located at the proximal hub of the balloon catheter and apply negative pressure for about 15 seconds until no bubbles appear in the contrast medium solution. Return to atmospheric pressure allowing contrast medium flow into the catheter lumen. Remove the syringe leaving a meniscus of contrast in the hub of the balloon lumen.
- 3 Remove air bubbles from theinflation device following the manufacturer's instructions.
- 4 Using stopcock attach inflation device to the PTA balloon catheter.

Avoid aire entering the system.
CAUTION: Avoid manipulation of the balloon catheter during removal from packaging and flushing of guidewire lumen. The GENOSS PTA Balloon Catheter is provided with a hydrophilic coating. Prior to inserting the catheter, activate the coating by immersing the catheter in normal saline, for approximately 30-60 seconds, or wiping down the catheter shaft with a saturated gauze sponge.
CAUTION: Do not wipe down the catheter surface with dry gauze.

H INSERTION TECHNIQUE

Access to the vessel must be sufficiently patent, sufficiently recanalized to proceed with further intervention. In case the guidewire has not already been positioned across the lesion, under fluoroscopy, insert a 0.035"(0.89mm) guidewire across the lesion following standard PTA techniques. Insert the proximal end of the guidewire into the distal tip of the catheter. Carefully insert the PTA balloon catheter through the introducer sheath and, under fluoroscopic guidance, carefully advance it through the peripheral vasculature following the guidewire towards the lesion.

CAUTION: Do not advance the guidewire or the balloon dilatation catheter if resistance is met without first determining the cause of resistance and taking remedial action. Position the balloon within the lesion using the two radiopaque markers located under the balloon as reference points for precise placement across the target lesion. To avoid kinking, advance the dilatation catheter slowly, in small increments until the proximal end of the guidewire emerges from the catheter.

I BALLOON INFLATION

CAUTION: Use only a mixture of sontrast medium and saline solution to fill the balloon(1:1).

Do not use any gaseous substances as a balloon inflation medium. Make sure all air bubbles are removed from the catheter and the inflation device before balloon inflation. Open stopcock on inflation device. inflate the dilation balloon gradually to dilate the lesion to the calculated diameter in accordance with the Compliance Chart. Apply a constant pressure for about 30 seconds.

The standard inflation medium is a 1:1 mixture of contrast medium and normal saline.

CAUTION: Do not exceed Rated Burst Pressure (RBP).

J BALLOON DEFLATION AND REMOVAL

Deflate the balloon catheter in accordance with standard PTA procedures. Apply negative pressures to the balloon for about 60-120 seconds. Ensure the balloon is fully deflated before carefully removing the catheter from the vessel.

K HOW SUPPLIED

The GENOSS PTA Balloon Catheter is supplied sterile and intended for single use only. The GENOSS PTA Balloon Catheter is sterilized by ethylene oxide gas. It will remain sterile as long as the packaging remains unopened and undamaged. Use product prior to labeled Expiration Date.

CAUTION: Do not use if the inner package is opened or damaged.

L STORAGE

Store at controlled room temperature (1°C~30°C) in a dry place.

M SYMBOLS

REF	Catalogue number	LOT	Batch code
	Do not re-use		Caution
	Do not use if package is damaged		Temperature limit
STERILE	Sterilized using ethylene oxide		Keep away from sunlight
	Use by date		Date of manufacture
	Consult instructions for use		Manufacturer
	Do not resterilize		Keep dry

N SIZE

Catheter length: 50cm, 80cm, 130cm

Diameter	Balloon length					
	20mm	40mm	60mm	80mm	100mm	120mm
4.0mm	•	•	•	•	•	•
5.0mm	•	•	•	•	•	•
6.0mm	•	•	•	•	•	•
7.0mm	•	•	•	•	•	•
8.0mm	•	•	•	•		
9.0mm	•	•	•	•		
10.0mm	•	•	•			
12.0mm	•	•	•			