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Bone substitute of GENOSS – **ORTHOPEdic OSTEON™ 3**

A Product Description

ORTHOPEdic OSTEON™ 3 is a synthetic osteoconductive bone graft substitute composed of 60% hydroxyapatite (HA) and 40% beta-tricalcium phosphate (β-TCP). ORTHOPEdic OSTEON™ 3 presents well interconnected porous structure, similar to that of human cancellous bone. It is supplied sterile by gamma irradiation. ORTHOPEdic OSTEON™ 3 is available as irregular shaped particle for orthopedic surgery.

B Indication for use

The use of ORTHOPEdic OSTEON™ 3 may be considered when autogenous bone is not indicated, or insufficient in quantity to fulfill the needs of the proposed surgical procedure. Case selection and use of appropriate surgical procedure are at the discretion of the physician.

- Bone filling
- Fractures with bone defects
- Pseudoarthrosis with or without bone defects
- Tibial osteotomy
- In certain cases or arthroplasty revision

C Contraindication

ORTHOPEdic OSTEON™ 3 should not be used in patients with:

- Chronic high dose therapy with corticosteroids
- Severe renal dysfunction, liver disease
- Established or imperfectly dried out osteomyelitis
- Necrosis at the recipient site
- Degenerative bone disease in the limb to be treated
- Intra-articular implantations
- Opening meniscus
- In the case of certain treatments under way, with corticoids or with drugs acting on the phosphocalcic metabolism, the use of ORTHOPEdic OSTEON™ 3 must be very carefully evaluated.

D Directions for use

- After exposure of the bony defect with mucoperiosteal flap, all inflammatory tissue and debris must be carefully removed using curettes.
- The surgical site should be cleaned with sterile water or sterile saline and aspirated to prepare the bone for the placement of ORTHOPEdic OSTEON™ 3.
- Mix ORTHOPEdic OSTEON™ 3 with sterile normal saline or with patient's blood in a sterilized dish.
- In order to enhance the formation of new bone, ORTHOPEdic OSTEON™ 3 should only be placed in direct contact with well vascularized bone. Cortical bone should be mechanically perforated.
- Transfer ORTHOPEdic OSTEON™ 3 granules to the defect site using a sterile instrument. Fill in the defect site loosely with ORTHOPEdic OSTEON™ 3 material to avoid crushing the particles and loss of trabecular architecture.
- Over filling of the defects should be avoided.
- The mucoperiosteal flaps should be sutured to achieve primary closure. A surgical dressing may be placed over the wound for one to two weeks.
- ORTHOPEdic OSTEON™ 3 is bone void filler without initial mechanical properties. Rigid fixation techniques may often be recommended.
- Site grafted with ORTHOPEdic OSTEON™ 3 should be allowed to heal approximately 6 months prior to implant placement.

E Precaution

- Effect on pediatric patients is not known.
- Effect on patients with a preexisting disease condition is not known.

F Warning

- Single use only. Do not resterilize or reuse.
- Not intended for immediate load-bearing.
- Do not overfill defects.
- Do not leave defect open.
- Do not compromise blood supply to the defect area.
- Do not apply excessive force when implanting ORTHOPEdic OSTEON™ 3 at the surgical site

- Do not use if package is opened or damaged or if expiration date has been exceeded.
- It is recommended to cover the implant site with a membrane when using ORTHOPEdic OSTEON™ 3.
- Use these products at the enough blood supply and bone contact site.
- Product should be protected from contamination due to touching patient's tongue or saliva.
- Open the cap from the arrow mark. (V)

G Adverse effect

No adverse reactions have been reported.

- Local complication could appear such as edema, numbness, bleeding after surgery, ulcer of soft tissue, infection etc.
- Numbness by injured nerve after surgery could appear, normally it is temporary phenomenon but might cause permanent damage.

H Symbols

	Catalogue number		Batch code
	Do not reuse		Date of manufacture
	Caution		Use by date
	Temperature limit		Manufacturer
	Sterilized using irradiation		Fragile, handle with care

I Storage

Store in room temperature away from sunlight(1~30°C)

J Expiration data

The expiry date is described on the container and packaging.

K Caution

Federal law restricts this device to sale by or on the order of a dentist or physician.

L Articles

REF	Size(mm)	Volume(cc)
O3G203005		5.0
O3G203010	2.0	10
O3G203015	~	15
O3G203020	3.0	20
O3G305005		5.0
O3G305010	3.0	10
O3G305015	~	15
O3G305020	5.0	20
O3G507005		5.0
O3G507010	5.0	10
O3G507015	~	15
O3G507020	7.0	20

Other sizes are available.

MS-P1691/IFU-O1302(Rev.2,2603)

GENOSS Co., Ltd.

Head Office: 12F, 76, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, Republic of Korea,
1201-1207, Gwanggyo Business Center, 156, Gwanggyo-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, Republic of Korea
Tel. +82-31-888-5100 www.genoss.com