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## 9 Electromagnetic Compatibility for TT-1.0

### Guidance and Manufacturer's Declaration – Electromagnetic Compatibility

Guidance and manufacturer's declaration – electromagnetic emissions		
TT-1.0 is intended for use in the electromagnetic environment specified below. The customer or the user of the TT-1.0 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-Guidance
RF emissions CISPR 11	Class 1	The TT-1.0 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group B	The TT-1.0 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded.  Warning: This equipment/system may cause radio interference of may disrupt the operation of nearby equipment. It may be necessary to take mitigation. Measures, such as re-orienting or relocating the TT-1.0 or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC 61000-3-3	Complies	

### Guidance and manufacturer's declaration – electromagnetic immunity

The TT-1.0 is intended for use in the electromagnetic environment specified below. The customer or the user of the TT-1.0 should assure that it is use in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humid it should be at least 30%.
Electrical fast transient/burst EN 61000-4-4	2 kV (for power supply lines) $\pm 1$ kV (for input/output lines)	$\pm 2$ kV (for power supply lines) $\pm 1$ kV (for input/output lines)	N/A
Surge EN 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV differential mode $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV common mode	$\pm 0.5$ kV, $\pm 1$ kV differential mode $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV common mode	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	0% $U_T$ (100% dip in $U_T$ ) for 0.5/1 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25/30 cycles <sup>a</sup> 0% $U_T$ (100% dip in $U_T$ ) for 250/300 cycles <sup>a</sup>	0% $U_T$ (100% dip in $U_T$ ) for 0.5/1 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25/30 cycles <sup>a</sup> 0% $U_T$ (100% dip in $U_T$ ) for 250/300 cycles <sup>a</sup>	N/A
Power frequency (50/60 Hz) Magnetic field EN 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE:  $U_T$  is the main voltage (AC) prior to the application of the test level.

<sup>a</sup> For example, 10/12 means 10 cycles at 50 Hz or 12 cycles at 60 Hz.

### Recommended separation distances between portable and mobile communication equipment and Derma LD Scalp Care

The TT-1.0 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TT-1.0 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TT-1.0 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter [m]				
	IEC 60601-1-2: 2007			IEC 60601-1-2: 2014	
0.01	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 2.7 GHz $d = 2.0\sqrt{P}$
	0.12	0.12	0.23	0.12	0.20
0.1	0.38	0.38	0.73	0.38	0.63
1	1.2	1.2	2.3	1.2	2.0
10	3.8	3.8	7.3	3.8	6.3
100	1.2	12	23	12	20

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

\*NOTE 1) At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

\*NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Guidance and manufacturer's declaration - electromagnetic immunity

The TT-1.0 is intended for use in the electromagnetic environment specified below.  
The customer or the user of the TT-1.0 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment- guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz – 80 MHz	3 Vrms	Portable and mobile RF communications equipment should not be used closer to any part of the Ultrasound System, including cables, than the recommended separation distance. This is calculated using the equation applicable to the frequency of the transmitter.  Recommended Separation Distance $d = 1.2\sqrt{P}$  IEC 60601-1-2:2007  $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5 GHz
	Outside ISM Bandsc amateur radio bands Bandsd  6 Vrms 150 kHz – 80 MHz In ISM bandsc amateur radio bands Bandsd		IEC 60601-1-2:2014 $d=2.0$ 80 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,a should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with following symbol:
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.7 GHz	6 V/m	

\* NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

\* NOTE 2: These guidelines may not apply in all situations, Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy
- b. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TT-1.0 is used exceeds the applicable RF compliance level above, the TT-1.0 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the TT-1.0. When the frequency range exceeds 150 kHz – 80 MHz, the electric field strength should be not higher than 3 V/m.
- c. The ISM (Industrial, Scientific and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz
- d. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

#### Guidance and Manufacturer`s Declaration – Electromagnetic Immunity

This product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipments should be used no closer than 30 cm (12 inches) to any parts of the product. Otherwise, the performance of this equipment could be degraded.

#### Guidance and manufacturer's declaration - electromagnetic immunity

The TT-1.0 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30cm (12 inches) to any part of the TT-1.0. Otherwise, degradation of the performance of this equipment could result.

Immunity test	Band a	Service a	Modulation	IEC60601 test level	Compliance level
Proximity fields from RF wireless Communications IEC61000-4-3	380 - 390 MHz	TETRA 400	Pulse modulation 18Hz	27 V/m	27 V/m
	430 – 470 MHz	GMRS 460 FRS 460	FM $\pm 5$ kHz deviation 1 kHz sine	28 V/m	28V/m
	704 – 787 MHz	LTE Band13, 17	Pulse modulation 217 Hz	9 V/m	9 V/m
	800 – 960 MHz	GSM800:900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28 V/m	28V/m

Immunity test	Band a	Service a	Modulation	IEC60601 test level	Compliance level
Proximity fields from RF wireless Communications IEC61000-4-3	1700 – 1990 MHz	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1,2,4,25 UMTS	Pulse modulation 217 Hz	28 V/m	28V/m
	2400 – 2570 MHz	Bluetooth WLAN 802.11b/g/n RFID 2450 LTE Band 7	Hz	28 V/m	28V/m
	5100 – 5800 MHz	WLAN 802.11a/n	Pulse modulation 217 Hz	9 V/m	9 V/m

\*NOTE : If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m. The 1m test distance is permitted by IEC 61000-4-3.

- a. For some services, only the uplink frequencies are included.
- b. The carrier shall be modulated using a 50% duty cycle square wave signal.
- c. As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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